

AirPex

Wireless Apex Locator

International Patent No. : PCT/CN2019/110429





Feel the Charm of Wireless Charging

Lay it down, Charge it up
Charging is as simple as placing
the AirPex on the charging base



True Wireless Apex Locator

Get rid of the wire to the fullest
Use without magnet wire
Accuracy through wirelessness



Model	AirPex
Dimensions	13cm x 11cm x 8cm ±1cm (Outer box)
Weight	0.35kg ±10%
Power supply	Lithium ion battery: 3.7V, 120mAh, ±10%
Charger power supply	AC 100-240 V, ±10%
Charger power output	5V --- 1A
Frequency	50/60Hz, ±10%
Power rating	<1W
Degree of protection	IPX0
Electrical safety class	Class II
Relative humidity	<80%
Ambient temperature	-20 °C ~ +55 °C

Air



Sky



Sunrise



AirPex
Wireless Apex Locator