



# Curing Pen


4 LEDs curing light






Smart controlling chip  
keeps light output consistent

Light intensity:  
 $1000\text{mW}/\text{cm}^2$ ,  $1500\text{mW}/\text{cm}^2$ ,  $2300\text{mW}/\text{cm}^2$   
 $600\text{mW}/\text{cm}^2$  (detect mode)



Wavelength: 385 - 515nm  
3 Blue, 1 Violet, 10w LED

360° rotatable head




Intensity digital display




Multiple curing options

### Ramp mode



### Pulse mode



### Caries detection



## Technical parameter:

Dimension: 24\*21\*230(mm)

Net weight: 109g

Light source head size:  $\phi 13\text{mm}$

Light output lens size:  $\phi 9.3\text{mm}$

Power supply: >1200 times

2300mW/cm<sup>2</sup>, 3 Sec

## Curing Pen

4 LEDs curing light

