

# Fast-Fill

Obturation Fill Handpiece



### Motorized GP Flow

- Provides consistent flow of GP at desirable rate
- Precludes voids and makes obturation homogeneous
- Eliminates hand fatigue that resulted from earlier gun type devices



### Fast-Fill Tips

- Autoclavable silver tips for multiple usages
- Silver tips maintain the thermoplastic temperature of GP and offers smooth flow inside canal
- 360° rotating tip for easy filling at all angles
- Extra long needle offers convenience and visibility
- Easy to maintain and replace



### Intuitive Design

- Sleek and ergonomic hand-piece
- Simple push button to dispense GP
- LED indicator reflects when the device is ready for use



### Comprehensive Display

- Three customizable memory programs - T1 to T3
- Temperature range 100°C - 200°C
- GP flow rate - Low, Medium, High
- Reserve GP level indicator

T1 200°C  
GP 80%

### Powerful Battery

- Highly efficient 2600mAh Lithium battery
- Direct Plug & Play - In case of low battery, use the device directly with adaptor



23 Gauge  
**Large**



25 Gauge  
**Small**



## Specifications

Heating Time	20 seconds to 200°C
Charging Time	4 hours to complete charge
Hexagonal Wrench	Use the wrench to bend the tip, rotate the tip and clean the device

*Fast-Fill*

