

Fast-Fill Fast-Pack pro



— Cordless Obturation System —

Fast-Fill



Useable cartridge tip,
economical design
GP pellet refill into device
Easy maintenance

Motorized gutta-percha
Push & stop, better flow
controlling

Digital displaying of rest
gutta percha

Real-time temperature
displaying
Set from 100~200°C

Powerful battery

3 ranges pushing speed, low、 mid、 high

Model	Fast-Fill
Dimensions	25cm x 10cm x 11cm±1cm (Outer box) Diameter:2.5mm-2.8mm Length: 14mm-16mm Applicable temperature: 100°C-200°C
Weight	1.1kg±10%
Power supply	Lithium ion battery: 3.7V, 2600mAh, ±10%
Charger power supply	AC 100-240 V, ±10%
Charger power output	6V---3A
Frequency	50/60Hz, ±10%
Temperature	100°C-200°C
Electrical safety class	Class II

| Obturation Fill Handpiece

Fast-Pack pro



4-5mm tip heating, safe and fast



Ring switch, 360°touchable

2 166°C

Real-time temperature displaying

Temperature 90°C

Customized temperature from 90~300°C

Model	Fast-Pack pro
Dimensions	20cm x 10cm x 11cm±1cm (Outer box)
Weight	1kg±10%
Power supply	Lithium ion battery: 3.7V, 2600mAh, ±10%
Charger power supply	AC 100-240 V, ±10%
Charger power output	6V---3A
Frequency	50/60Hz, ±10%
Power rating	15VA
Temperature	90°C-300°C
Electrical safety class	Class II

| Obturation Pack Handpiece

Fast-Fill Fast-Pack pro

