

HyperLight-G

www.eighteeth.com



HYPERLIGHT-G



Receptor

Use with film, digital sensor, phosphor plate




Large Battery Capacity

More than 300 exposures in one full charge

Preset Parameter

Divided into adults and children, and can be set and controlled for different teeth



Smart angulation system, faster your workflow

Locking system





Safe

- Internal x-ray shielding
- Locking system to avoid unintentional exposure
- Backscatter shield for operator and patient protection(*optional)

Portable

- Cordless, handheld and extremely mobile, move freely between rooms, locations and cities

Economical

- Use with film, digital sensor or phosphor plate
- Prolonged working life, at least 80,000 exposures

HyperLight-G

www.eighteeth.com



| | |
|---------------------|------------------------|
| Model | : HyperLight-G |
| Tube Voltage | : 65 kV |
| Tube Current | : 2 mA |
| Exposure Time Range | : 0.02-2 Sec |
| Weight | : 1.9 kg |
| Battery | : 2000mAh X4 |
| Line Voltage | : 100V-240V, 50Hz/60Hz |